


UV752/ UV2000 Spectrophotometer



| Technique parameter | |
|------------------------------------|--|
| Optical System | Single beam, C-T type, 1200 lines/mm grating |
| Detector | Silicon photo diode |
| Monitor | 4.3 inch touch screen display |
| Light Source | Tungsten lamp (20W / 12V) Deuterium lamp In the 340 nm automatic switching |
| Spectral Bandwidth | 4nm |
| Wavelength Range | 200~1000nm (step interval 0.1nm) |
| Wavelength Accuracy | ±2.0nm (Auto-calibration when starting) |
| Wavelength Repeatability | ≤1.0nm |
| Stray Light | ≤0.5%T (at 220,340nm) |
| Photometric Range | T: -1.0~200.0% , A: -0.5~3.000Abs, F: 0~9999, C: 0~9999 |
| Transmittance Accuracy | ±0.5%T (0~100%T) |
| Transmittance Repeatability | ≤0.2%T |
| Stability | Dark current < 0.2%T , Bright current < 0.5%T |

Address: 5-6/F, Building 8, 456 Hong Cao Road, Xu Hui District, Shanghai 200233 P. R. China

Tel: 86-21-64951509

E-mail : info@hengping.com

Sales Hotline: 86-21-64951010

E-mail : sales@hengping.com

Customer Service: 86-21-64951509

<http://www.hengping.com>